



Symmetry Dental Equipment
Symmetry 4.0 Collection
SYMBDS E - Rear Mount Floating Worksurface with Duo Delivery Package
Specification Page

Integrates with all Symmetry 4.0 Rear Treatment Consoles



Decor Finish for Door

Cabinet Sold Separately

Dr. Set Up:

- Three HP automatic
- One 3 way air / water syringe
- Wet / Dry foot control
- Silcryn tubing
- 1/4" water quick disconnect
- Easy to read pressure guage
- Pressure adjustment for each handpiece
- Non-retracting water coolant, adjustmetnt, and purge
- Master on / off toggle with indicator
- Master valve assembly with filter and regulator
- 7' umbilical tubing

Assistant Set Up:

- Dual autoclavable HVE
- Saliva Ejector
- One 3 way air / water syringe
- Silcryn tubing
- Vacuum canister with removable screen

Integration Set Up:

- Matching solid surface worksurface for Symmetry 4.0 series
- Cabinet mounted water bottle system

Decor Finish for Door	
Counter Finish	